



# High Frequency OSCillation in ARDS

Dr Elena Lambert

Clinical Fellow for the OSCAR Trial

# Introduction

- What is HFOV
- About the OSCAR Trial
- The Vision Alpha ventilator
- Day to day use of the ventilator
- Protocols
- Support
- Summary and Questions

# Agenda



- Theory
  - This morning
- Practical
  - This afternoon
  - Setting up the ventilator
  - SimMan demonstration
- Q&A
- Feedback

# What is HFOV?

- 'CPAP with a wiggle'
- High frequency ventilation
- Tidal volumes approach or are less than anatomical dead space
- Lungs kept 'open' with a constant airway pressure
- Allows clinician to deal with oxygenation and ventilation separately

# History of HFOV

- Initial concept, 1915, smoke blown down a tube travels in a very thin column, could you ventilate with only a small volume?
- Trialled in animals from 1975 but results not replicated in human model
- Used successfully in neonates and then Paediatrics since 1983
- Several small trials in adults since 1999
- Used in several units in UK but often 'when all else fails'

# Theory behind HFOV

\*Weavind and Wenker 2000



- Bulk flow can still provide conventional gas delivery to proximal alveoli with low regional dead space volumes.
- Coaxial flow. Gas in the centre flows inward, while gas on the periphery flows outward. This can develop because of the asymmetric low profile of high velocity gases.

# Theory behind HFOV (2)

\*Weavind and Wenker 2000



- Pendelluft can mix gases between lung regions having different impedances.
- Taylor dispersion can produce a mixing of fresh and residual gas along the front of a flow of gas through a tube.

# Theory behind HFOV (3)


\*Weavind and Wenker 2000







- Augmented molecular diffusion can occur at the alveolar level secondary to the added kinetic energy from the oscillations.

# Why HFOV?

- We know that 'conventional ventilation' can be harmful
- ARDS has high morbidity and mortality rates
- Lung protective ventilation
- Theorised to reduce biotrauma





# The OSCAR Trial

- High frequency OSCillation in ARDS
- Several small studies comparing CV and HFOV – but not enough evidence
- Need to know if it is better
- 10 centre, open, randomised control trial
- Early ARDS rather than as a rescue therapy
- Start recruiting 1/11/07, for 3 years
- Target = 2 patients per month per unit

# Other Trials

- Two other current trials
- USA and Canada
- Similar criteria to ours
- Small
- Canada has view to extend if preliminary results are good
- Using Sensormedics ventilator

# Inclusion Criteria



- $\geq 16$  years
- $\geq 35$  kg
- ARDS
  - $\text{PaO}_2/\text{FiO}_2$  ratio  $\leq 26.7$  kPa, with a PEEP  $\geq 5$  cmH<sub>2</sub>O on 2 ABG's, 12 hours apart
  - Bilateral infiltrates on CXR
  - 1 or more RF for ARDS
- Receiving artificial ventilation
- Likely to require  $\geq 48$  hours of ventilation
- 'Substantially uncertain' as to the utility of HFOV for this patient
- Ventilated for  $< 7$  consecutive days (168 hours)

# Exclusion Criteria



- LA hypertension
- Moderate - severe airway disease causing airflow limitation
- Lung biopsy or resection (this admission)
- Patient refuses
- Patient unable to consent (language)
- If active treatment is/to be withdrawn
- In another trial (30 days prior to randomisation)
- Previously enrolled in the OSCAR trial
- Any other condition likely to make HFOV hazardous

# Vision Alpha Ventilator



- Both HFOV and CV
- New to this country
- Made in Japan
- Currently used in Germany
- CE marked
- Motorbike muffler to reduce noise

# Front view




- Divided into main (upper) and blower (bottom) units
- Operation panel
- Respiratory circuit connector sections
- HFOV drive section
- Dedicated humidifier bracket

# Rear View




- Interface block, currently only used for service
- Gas connector section
- Breaker section
- Filter section
- Auxiliary breaker section
- No battery

# Patient set



- Humidifier
- Expiratory filter
- Circuit
- Oscillator diaphragm
- Disposable

# Circuit




- No rain out
- Narrower tubing
- Heated from the outside
- Disposable

# Impedance valve



- Expensive!
- Requires EOG sterilisation between patients

# Expiratory Filter



- Heated so does not require changing as frequently
- Easy to change

# Humidifier



- Opaque
- Large membrane
- Auto fills
- 30 day life

# Control Panel



- A: CV controls
- B: FiO<sub>2</sub> control
- C: Mode adjustment
- D: Accept button
- E: HFO controls
- F: Screen
- G: LED bar graph
- H: Silence alarm

# HFOV Controls



# Conventional Ventilation Controls




# FiO<sub>2</sub> and Mode Controls




# Accept, Silence and LED Bar




# HFO Screen




# CV screen

Current mode: P C V-AC		Normal operation	Alarm sound ON												
Inspiratory pressure upper limit	35	cmH <sub>2</sub> O													
Inspiratory pressure lower limit	3	cmH <sub>2</sub> O													
PEEP lower limit	0	cmH <sub>2</sub> O													
Mandatory ventilation volume lower limit	0	mL													
Spontaneous ventilation volume lower limit	0	mL													
Breath rate upper limit	150	SPM													
Minute-volume-of-ventilation upper limit	60.0	L													
Minute-volume-of-ventilation lower limit	1.00	L													
Apnea monitoring time	20	sec													
<table border="1"> <thead> <tr> <th colspan="2">IMV patient data</th> </tr> </thead> <tbody> <tr> <td>Peak inspiratory pressure</td> <td>9.1 cmH<sub>2</sub>O</td> </tr> <tr> <td>Mean airway pressure</td> <td>4.3 cmH<sub>2</sub>O</td> </tr> <tr> <td>Breath rate</td> <td>20 B<sub>P</sub>M</td> </tr> <tr> <td>Minute volume of ventilation</td> <td>3.3 L<sub>M</sub></td> </tr> <tr> <td>Tidal volume</td> <td>166 mL</td> </tr> </tbody> </table>				IMV patient data		Peak inspiratory pressure	9.1 cmH <sub>2</sub> O	Mean airway pressure	4.3 cmH <sub>2</sub> O	Breath rate	20 B <sub>P</sub> M	Minute volume of ventilation	3.3 L <sub>M</sub>	Tidal volume	166 mL
IMV patient data															
Peak inspiratory pressure	9.1 cmH <sub>2</sub> O														
Mean airway pressure	4.3 cmH <sub>2</sub> O														
Breath rate	20 B <sub>P</sub> M														
Minute volume of ventilation	3.3 L <sub>M</sub>														
Tidal volume	166 mL														
IMV setup	IMV alarm	IMV data	Waveform												
Loop		H F O													

# Waveform screen



# Loop screen



# Differences between CV and HFOV

- Tidal volumes  $\leq$  dead space
- Decreasing the frequency (i.e. RR) helps remove CO<sub>2</sub>
- Different values to observe and record
- Auscultation of little value
- No rain-out in tubing owing to new humidifier

# Family

- Possibility of more concerns as patient looks less 'alive' on HFOV compared to CV
- Or may think that this is a 'miracle machine' that will save patient
- Family need to be informed and not given unrealistic hope

# Daily data collection




- Mean Airway Pressure
  - A control of oxygenation
  - Start at value according to protocol
  - Increase to promote recruitment and oxygenation
  - Can be decreased when happy with oxygenation to a more constant level
  - Value on screen will fluctuate slightly

# Daily data collection




- Frequency ( $H_z$ )
  - A control of ventilation
  - Start at value according to protocol
  - Decrease this to remove  $CO_2$
- Cycle volume (ml)
  - A control of ventilation
  - Start at value according to protocol
  - This will change with a change in frequency
  - Increase this to remove  $CO_2$
  - Affects the Amplitude

# Daily data collection



- Measured Amplitude (cmH<sub>2</sub>O)
  - Reflects the compliance of the circuit/lungs
  - Changes can be subtle but can be effective warning system
  - Will change if the frequency is altered
  - Unexplained rise: Secretions/Obstruction
  - Unexplained fall: Pneumothorax

# Daily Data Collection



- Base (Bias) Flow
  - Usually stable
  - Set at level according to protocol
  - May need to be increased e.g. if patient has lighter sedation

# Suctioning

- Always suction patient before going on to HFOV
- Suctioning possible with closed system device
- Need to be aware that this can cause de-recruitment
- Low pressure alarm may sound, PaO<sub>2</sub> may drop
- May need to re-recruit lungs after as algorithm
- Physio: possible as normal

# Problems at initiation

- Coughing
- MAP swings
- Alarms
- Hypotension – most common cause at this stage is hypovolaemia


# Problems once established

- Pneumothorax
  - Hypotension (acute), fall in PaO<sub>2</sub>, asymmetric chest wall movement
- Air trapping
  - Hypotension (gradual), fall in PaO<sub>2</sub>, rise in CVP
- Blocked ET Tube
  - Fall in PaO<sub>2</sub>, rise in PaCO<sub>2</sub>, increase in Amplitude

# HFOV Protocol


## Control of Oxygenation

- Step 1 = Recruitment
- Step 2 =  $FiO_2$  delivery
- Step 3 = Reducing  $P_{AW}$



# HFOV Protocol

## Management of Ventilation



© Iain Mackenzie, Bill Tunnicliffe, Sanjoy Shah and the OSCAR trial management group, 2007.

# CV Protocol

- As ARDS 'best practice'
- Patients ventilated using tidal volumes of 6-8ml kg<sup>-1</sup> body weight
- Plateau pressure of no greater than 30cmH<sub>2</sub>O.
- The ventilation target = normal arterial pH
- Permissive hypercapnia to avoid plateau pressures above 30cmH<sub>2</sub>O.
- Oxygenation target = PaO<sub>2</sub> of 7.5-11.5kPa or an SpO<sub>2</sub> of 88–95%.

# CV Protocol

- The recommended combinations of  $\text{FiO}_2$  and PEEP are:

$\text{FiO}_2$	PEEP
0.3	5
0.4	5
0.4	8
0.5	8
0.5	10
0.6	10
0.7	10-14
0.8	12-14
0.9	12-16
1.0	12-18

# Support

- Me
- Management Group
- Inspiration Healthcare/Novalung
- Website
- Abbreviated manual attached to ventilator

# Summary


- HFOV – frequency 5-15Hz, tidal volumes  $\leq$  dead space, lungs kept 'open'
- Separates oxygenation and ventilation
- Protocols guide management of patient
- Several of sources of support
- Recruiting from 1<sup>st</sup> November 2007
- 2 patients per centre per month (min)

# Practical

- How to assemble
- SimMan Demonstration

# Any Questions?

(Feedback Forms)



# My Contact Details

Dr Elena Lambert

Clinical Fellow OSCAR Trial

ICS Trials Group

Kadoorie Centre

John Radcliffe Hospital

Oxford. OX3 9DU

- Tel: 01865 221755
- Fax: 01865 857611
- Work Mobile:  
07825733953
- Bleep: 1114 (via John Radcliffe switchboard)
- Email:  
[OSCAR.Trial@nda.ox.ac.uk](mailto:OSCAR.Trial@nda.ox.ac.uk)